



ease@T

SOFT, SAFE, SYNCED BREAST HEALTH

## CO REVIEWER USER MANUAL



**Login to your account**  
Enter your email below to login to your account

Email



Password

[Forgot your password?](#)

**Login**

Don't have an account? [Sign Up](#)

They can sign in using the email and password provided by the scan centre.




**Verification Code**  
We have sent the verification code to  
**sperformingprovider@gmail.com**

00:57  
Waiting for Passcode...

**Verify**


A verification code will be generated and sent for verification.



The image shows a 'User Consent Form' for the 'EaseQT Platform'. The form is presented in a window with a pink header. It contains four sections of consent text:

- 1. Consent to Platform Use**  
I hereby acknowledge that I am voluntarily registering to use the EaseQT platform for purposes directly related to clinical imaging, workflow management, or reporting.
- 2. Data Privacy & Security**  
I understand that my personal information and any patient-related data accessed through this platform may be securely stored and processed in compliance with global privacy laws (including HIPAA, DPA 2021, and GDPR where applicable). I agree not to share, download, or misuse any protected information obtained via the platform.
- 3. Confidentiality Agreement**  
I agree to maintain full confidentiality of any data I access, including but not limited to patient records, reports, clinical communications, and platform operations. Any breach of confidentiality may result in legal and professional consequences.
- 4. Data Consent**  
I consent to the transfer, storage, and processing of data as necessary for platform function and continuity of care. I understand that all such data will be handled securely and only by authorized personnel.

The consent form must be filled by the co reviewer.




The image shows a 'Reset Your Password' form. On the left is a collage of medical-related images. On the right, the form has a pink background and includes the following elements:

- easeQT** logo with the tagline 'SOFT. SAFE. SYNCED BREAST HEALTH'.
- Reset Your Password** title.
- E-Mail: `sperformingprovider@gmail.com`
- New Password** input field with a strength indicator.
- Password requirements (all checked):
  - At least 8 characters
  - Contains a number
  - Special character (!@#%&^&\*)
- Confirm Password** input field with a strength indicator.
- Reset Password** button.


The password must be reset for the security purpose.

## DASHBOARD



In dashboard, the scan centre co reviewer can view patient brochure, patient consent, tech guideline and tech consent.

## PATIENT QUEUE



In the patient queue, the scan centre performing provider can view the patient list, patient forms, technician view and DICOM files. They can fill out reports, check the report status, assign the case to another user and add pending remarks if needed and can view previous remark also. Click on **View** to proceed with filling out the report. And they can multi select and download DICOM and Final reports. They can upload and download patient's old report.

Patient Form

Technician Form

Dicom

Report

Final Report

No Reports Found

Other Old reports

S.No	Report Title	File
No Reports Found		

Time stamp

Created by	Start	End	Status
Super Admin	19 Sept 2025 20:55	19 Sept 2025 20:56	Technologist Form Fill
Super Admin	19 Sept 2025 20:56	19 Sept 2025 20:57	Reviewed 1
MSC Co Reporting Doctor	19 Sept 2025 20:57	19 Sept 2025 21:21	Reviewed 2
MSC Co Reporting Doctor	19 Sept 2025 21:21	-	-

Status saver

Load Template

Reset to Default

Predraft

Draft

Reviewed 1 Correct

Reviewed 1 Edit

Reviewed 2 Correct

Reviewed 2 Edit

Insert Signature

Sign Off

Normal

NAME	U S	DOB	September 5, 1972 (53 y)
GENDER	F	SCAN CENTER	MSC, Bangalore
USERID	U	DATE OF VISIT	Sep 19, 2025

**QT ULTRASOUND BREAST IMAGING**  
 Indication: Screening - Routine Annual Checkup.  
 History:  
 Reason for having this QT scan (Patient Stated): 333.  
 Patient is asymptomatic.

**TECHNIQUE:** Transmission and reflection multipolar 3-dimensional ultrasound imaging of both breasts was performed using the QT Ultrasound Series 2000 Model-A scanner. Breast density was determined using the Quantitative Breast Density calculator. Images were reviewed in the QTViewer v2.0.2. The nipple-areolar complex, skin, Cooper's ligaments, breast fat distribution, penetrating arteries and veins, glandular and ductal tissues were evaluated. Images were reviewed in coronal, transaxial and sagittal planes.

**SYMMETRY:** symmetry  
 The breast examination indicates symmetrical size

Copyright © Wellthreen

They have access to Reviewed 2 correct and Reviewed 2 edit after reviewing the report. If co reviewer doctor needs to edit anything they can edit in final report.

Patient Form

Technician Form

Dicom

Report

Final Report

No Reports Found

Other Old reports

S.No	Report Title	File
No Reports Found		

Time stamp

Created by	Start	End	Status
Super Admin	19 Sept 2025 20:55	19 Sept 2025 20:56	Technologist Form Fill
Super Admin	19 Sept 2025 20:56	19 Sept 2025 20:57	Reviewed 1
MSC Co Reporting Doctor	19 Sept 2025 20:57	19 Sept 2025 21:21	Reviewed 2
MSC Co Reporting Doctor	19 Sept 2025 21:21	-	-

Status saver

Load Template

Reset to Default

Predraft

Draft

Reviewed 1 Correct

Reviewed 1 Edit

Reviewed 2 Correct

Reviewed 2 Edit

Insert Signature

Sign Off

Normal

NAME	U S	DOB	September 5, 1972 (53 y)
GENDER	F	SCAN CENTER	MSC, Bangalore
USERID	U	DATE OF VISIT	Sep 19, 2025

**QT ULTRASOUND BREAST IMAGING**  
 Indication: Screening - Routine Annual Checkup.  
 History:  
 Reason for having this QT scan (Patient Stated): 333.  
 Patient is asymptomatic.

**TECHNIQUE:** Transmission and reflection multipolar 3-dimensional ultrasound imaging of both breasts was performed using the QT Ultrasound Series 2000 Model-A scanner. Breast density was determined using the Quantitative Breast Density calculator. Images were reviewed in the QTViewer v2.0.2. The nipple-areolar complex, skin, Cooper's ligaments, breast fat distribution, penetrating arteries and veins, glandular and ductal tissues were evaluated. Images were reviewed in coronal, transaxial and sagittal planes.

**SYMMETRY:** symmetry  
 The breast examination indicates symmetrical size

Copyright © Wellthreen

There is an option called "Reset to Default." If the report was done incorrectly, they can use this option to revert it to the default state.

Patient Form

Technician Form

Dicom

Report

General

Right

Left

Impression + Reco

Final Report

← Back

Patient Name: sofia

Patient ID: P100016

Gender: Female

DOB: August 4, 1997 (28 Y)

Patient Form: S Form

Date: 18 Sept 2025

Scan Center: MSC

Other reports

S.No	Report Title	File
No Reports Found		

Other Old reports

S.No	Report Title	File
1	Thermogram	<a href="#">View</a>

Time stamp

Created by	Start	End	Status
sofia	18 Sept 2025	18 Sept 2025	Patient Intake For...

Final Report Preview

Ease QT 10.10 Auto Report ☒

Normal

B

I

U

NAME	sofia	DOB	August 4, 1997 (28 y)
GENDER	F	SCAN CENTER	MSC, Bangalore
USERID	P100016	DATE OF VISIT	Sep 18, 2025

QT ULTRASOUND BREAST IMAGING

Indication: Screening - Routine Annual Checkup.

History:

Reason for having this QT scan (Patient Stated): sss.

Patient is asymptomatic.

Imaging:

- sss.

TECHNIQUE: Transmission and reflection multiplanar 3-dimensional ultrasound imaging of both breasts was performed using the QT Ultrasound Series 2000 Model-A scanner.

Copyright © Wellthreen

Back

1/1

Next

z0SuYYZcbHUgC4gse+1qqr7Y0J7gUAEhCYITBAJM...

1 / 15

100%


1

S.No	Title	PAGE No
1	Introduction	2

Copyright © Wellthreen

The users can see uploaded Patient's old reports and their current report.

## ANALYTICS



In analytic, the co-reviewer can view their total case count, TAT, previous 6 months count, Reco code and artifacts.

## **GENERAL TIPS**

1. If you open something in edit mode, save it before exiting-otherwise, others won't be able to use it.
2. After addressing pending remarks, leave another remark saying "Done" and set Assign to as None.
3. The Filter option is available for sorting.
4. For bulk downloads, use the checkboxes to select patients and then download in bulk.
5. To add a pending remark, click the + icon.
6. When viewing a report, always click Back to exit-otherwise, the case you opened may freeze.
7. Use Image Preview to add images to a report.
8. To refresh the page, press Ctrl + Shift + R.
9. For DICOM downloads and voice dictation, use Chrome or Edge. (For Firefox, click Save. Do not use Internet Explorer.)
10. Always log out properly-do not just close the browser or tab.

## **Reviewer**

1. Review 2 is not reversible - admin intervention is required to undo it.
2. easeQT is disabled after signing off in R2.