


ease @T


SOFT, SAFE, SYNCED BREAST HEALTH

TECHNICIAN LOGIN USER MANUAL



Once the patient has completed the form, the technician can begin filling out the technician form.

The technician signs in using the email and password provided by Scan centre.



A verification code will be generated and sent for verification.



The image shows a 'User Consent Form' for the 'EaseQT Platform'. The form is presented in a light pink window with a close button in the top right corner. It contains four numbered sections: 1. Consent to Platform Use, 2. Data Privacy & Security, 3. Confidentiality Agreement, and 4. Data Consent. Each section includes a paragraph of text explaining the user's agreement to the platform's terms and conditions. The background of the window shows a collage of medical images, including ultrasound scans and a person in a clinical setting.

easeQT
SOFT, SAFE, SYNCED BREAST HEALTH

User Consent Form

EaseQT Platform


1. Consent to Platform Use
I hereby acknowledge that I am voluntarily registering to use the EaseQT platform for purposes directly related to clinical imaging, workflow management, or reporting.

2. Data Privacy & Security
I understand that my personal information and any patient-related data accessed through this platform may be securely stored and processed in compliance with global privacy laws (including HIPAA, DPA 2021, and GDPR where applicable). I agree not to share, download, or misuse any protected information obtained via the platform.

3. Confidentiality Agreement
I agree to maintain full confidentiality of any data I access, including but not limited to patient records, reports, clinical communications, and platform operations. Any breach of confidentiality may result in legal and professional consequences.

4. Data Consent
I consent to the transfer, storage, and processing of data as necessary for platform function and continuity of care. I understand that all such data will be handled securely and only by authorized personnel.

The technician will sign the consent form.



The image shows a 'Reset Your Password' form for the 'EaseQT Platform'. The form is presented in a light pink window. It features the EaseQT logo at the top, followed by the title 'Reset Your Password' and the email address 'E-Mail: techniciancentre@gmail.com'. Below this, there are two password input fields: 'New Password' and 'Confirm Password'. The 'New Password' field has a dropdown menu with options: 'At least 8 characters', 'Contains a number', and 'Special character (!@#%&*^&*)'. The 'Confirm Password' field has a dropdown menu with the option 'Confirm password'. At the bottom of the form is a 'Reset Password' button. The background of the window shows a collage of medical images, including ultrasound scans and a person in a clinical setting.

easeQT
SOFT, SAFE, SYNCED BREAST HEALTH

Reset Your Password

E-Mail: techniciancentre@gmail.com

New Password

Password

At least 8 characters
Contains a number
Special character (!@#%&*^&*)


Confirm Password

Confirm password

Reset Password


The technician needs to reset their password before logging in.

DASHBOARD



In the dashboard, the technician can view the patient brochure, technician guidelines, and consent form.

PATIENT QUEUE



In the patient queue, the patient list and patient form will be displayed. The technician needs to fill out the technician form and upload the DICOM files.

The screenshot shows the 'Dicom' upload section of the easeQT interface. It features two parallel upload zones. Each zone includes a header 'Upload Files' with a sub-instruction 'Select and upload the files of your choice'. Below this is a large dashed box with an upload icon and the text 'Choose a file or drag & drop it here ZIP files', accompanied by a 'Browse File' button. At the bottom of each zone, a file upload status is shown: 'Radiology (60).zip' (100487 KB of 100487 KB) is 'Uploading...' with a blue progress bar, while 'Radiology (64).zip' (4755 KB of 4755 KB) is 'Completed' with a green progress bar. A 'Submit' button is located at the bottom center of the upload area.


Once the technician completes the form, they will upload the patient's DICOM files.

The screenshot displays the 'Patient Queue' table. The table has columns for Date, Scan Centre, Patient Name, Form, Tech Form, DICOM (L, R), View, Old Report, QT Report, Report Status, Assign, Pending Remarks, and Remarks. The row for 'sofia' on 'Sep 18, 2025' is highlighted. The 'DICOM' column for this row shows '(1) ↓' and a red box highlights the 'View' column. The 'Report Status' for this row is 'DRAFT'. The table also includes a 'Download Dicom' button, a 'Download Final Reports' button, a 'Clear All Filters' button, and a search bar. The footer shows 'Rows per page 50', navigation buttons, and 'Page 1 of 1'.

	Date	Scan Centre	Patient Name	Form	Tech Form	DICOM L	DICOM R	View	Old Report	QT Report	Report Status	Assign	Pending Remarks	Remarks
<input type="checkbox"/>	Sep 17, 2025	MSC	TestingT2	Dc - View	View	(1) ↓	-	View DICOM		View	SIGNED OFF (A)	Assign	No remarks...	+
<input type="checkbox"/>	Sep 18, 2025	MSC	TestingT2	S - View	View	(1) ↓	-	View DICOM		View	REVIEWED 2	Assign	R1 R2 stuck	+
<input type="checkbox"/>	Sep 22, 2025	MSC	TestingT2	S - View	View	(1) ↓	-	View DICOM		View	SIGNED OFF (A)	Assign	No remarks...	+
<input type="checkbox"/>	Sep 18, 2025	MSC	sofia	S - View	View	(1) ↓	-	View DICOM		-	DRAFT	Assign	No remarks...	+

After uploading DICOM, the technician can view the uploaded DICOM files and can also download them if needed. The technician can upload multiple DICOM files at once. He technician can select multiple DICOM files along with the final report and download them all at once.

ANALYTICS



In Analytics, the technician can view their case count, TAT (Turnaround Time), scan indications, Artifacts and the case count for the previous six months. They can click on month to filter the data.

GENERAL TIPS

1. If you open something in edit mode, save it before exiting-otherwise, others won't be able to use it.
2. After addressing pending remarks, leave another remark saying "Done" and set Assign to as None.
3. The Filter option is available for sorting.
4. For bulk downloads, use the checkboxes to select patients and then download in bulk.
5. To add a pending remark, click the + icon.
6. When viewing a report, always click Back to exit-otherwise, the case you opened may freeze.
7. Use Image Preview to add images to a report.
8. To refresh the page, press Ctrl + Shift + R.
9. For DICOM downloads and voice dictation, use Chrome or Edge. (For Firefox, click Save. Do not use Internet Explorer.)
10. Always log out properly-do not just close the browser or tab.

Technician

1. Tech Form Rules: Once submitted, tech forms cannot be changed. Add any corrections or notes in Remarks instead.