



ease @T

SOFT, SAFE, SYNCED BREAST HEALTH

SCAN CENTRE ADMIN LOGIN USER MANUAL



easeQT
SOFT, SAFE, SYNCED BREAST HEALTH

Login to your account
Enter your email below to login to your account

Email
scanadmin@gmail.com



Password

[Forgot your password?](#)

Login

Don't have an account? [Sign Up](#)

The scan centre admin can sign in using email and password provided by Wellthgreen.



easeQT
SOFT, SAFE, SYNCED BREAST HEALTH

Verification Code
We have sent the verification code to
scanadmin@gmail.com

□ □ □ □ □

00:58
Waiting for Passcode...

Verify

A verification code will be sent to the email for verification and it needs to be entered in the application.

I consent to the transfer, storage, and processing of data as necessary for platform function and continuity of care. I understand that all such data will be handled securely and only by authorized personnel.

5. Credential Use and Responsibility

I confirm that all credentials, licenses, or certifications I provide during registration are accurate and current. I agree to keep my login credentials confidential and to notify the platform administrator of any suspected breach or unauthorized access.

6. Communication & Updates

I consent to receiving essential updates, notifications, and communications related to the EaseQT platform via email or in-app alerts.


☒ I agree to the above universal terms and policies.


7. Scan Centre Manager Consent

- I am responsible for managing users, cases, or workflows within my assigned center or group.
- I agree to uphold security and privacy protocols in all operations.
- I will not access or modify clinical data outside of my administrative scope.

☒ I agree to the Admin role-specific terms.

The consent form needs to be signed to proceed.



 **easeQT**
SOFT. SAFE. SYNCED BREAST HEALTH

Reset Your Password


E-Mail: scanadmin@gmail.com

New Password

At least 8 characters
Contains a number
Special character (!@#%&^&*)


Confirm Password

The scan centre admin needs to rest the password for security purpose.



In the dashboard, the scan centre can view the patient brochure, patient consent form, technician guidelines and technician consent form.

ADMINISTRATION



In Administration, the scan centre admin can add their technicians, performing providers, co-reporting doctors and Patients.

Dashboard
Administration
Patient Queue
Analytics
MSC Manager
CMSCM2
Scan Center Manager

← SCAN CENTRE

SCAN CENTER DETAILS

NAME : MSC

ADDRESS : Bangalore

TELEPHONE : 0433-3333333

WEBSITE : www.centre@com

EMAIL : testingmsc@gmail.com

SCAN CENTER CODE : MSC

Manage Technician
Manage Performing Provider
Manage Reviewer
MANAGE PATIENTS
Save

Copyright © Wellthgreen

Dashboard
Administration
Patient Queue
Analytics
MSC Manager
CMSCM2
Scan Center Manager

← MANAGE TECHNICIAN

+ Add Technician

Name	Technician ID	Mail ID	Mobile Number	Status	Manage Profile
MSC Technician	CMSCM3	msctechnician@gmail.com	9842653428	Active	Edit
Technician	MSCT2	technician@gmail.com	9876543219	Active	Edit
Rohith	MSCT1	rohithganesan208@gmail.com	9876543213	Active	Edit

Rows per page
5
<< < Previous 1 Next > >>
Page 1 of 1

Copyright © Wellthgreen

The centre manager can add and manage all users associated with their scan centre. The admin can add technician by click on ‘Add technician’.

Dashboard
Administration
Patient Queue
Analytics
scan admin CMSCM2

ADD TECHNICIAN

Add Photo

Personal Details

Full Name *
Enter Full Name

E-Mail *
m@example.com

Social Security Number *
Enter Social Security Number

Contact Number *
USA (+1) Enter Phone Number

Date Of Birth *
Select a date

Driver's License *
Choose File No file chosen

Professional Details

Scan Center Name *
MSC Gokul Scan Center

Digital Signature *
Choose File No file chosen

Upload License *
Choose Files No file chosen

☒ Trained in QT protocol

[< Back](#)
[Submit >](#)

The centre admin can add the personal and professional details required to add a technician. They can also attach the person's driver's license, and social security number.

Dashboard
Administration
Patient Queue
Analytics
MSC Manager CMSCM2 Scan Center Manager

MANAGE TECHNICIAN

[+ Add Technician](#)

Name	Technician ID	Mail ID	Mobile Number	Status	Manage Profile
MSC Technician	CMSCM3	msctechnician@gmail.com	9842653428	Active	Edit
Technician	MSCT2	technician@gmail.com	9876543219	Active	Edit
Rohith	MSCT1	rohithganesan208@gmail.com	9876543213	Active	Edit

Rows per page 5
<< < Previous 1 Next > >>
Page 1 of 1

Edit Technician

Technician ID: CMSCT3

Status: Active

Personal Details

Full Name: MSC Technician

E-Mail: msctechnician@gmail.com

Social Security Number: 1234567890

Contact Number: USA (+1) 9842653428

Date Of Birth: September 19th, 2002

Driving License: Upload Driving License

Driving License Document

Professional Details

Copyright © Wellthgreen

If the centre admin wants to change the information in a particular column, they can click the “Edit” option, make the necessary changes and submit. If they want to deactivate a user, they can change the user’s status from “Active” to “Inactive.”

The same process applies for scan centre performing provider, reviewer and technician.

ADDING PATIENTS

SCAN CENTER DETAILS

NAME: MSC

ADDRESS: Bangalore

TELEPHONE: 0433-3333333

WEBSITE: www.centre@com

EMAIL: testingmsc@gmail.com

SCAN CENTER CODE: MSC

Manage Technician Manage Performing Provider Manage Reviewer **MANAGE PATIENTS** Save

Copyright © Wellthgreen

Dashboard

Administration

Patient Queue

Analytics

M

MSC Manager

CMSCM2

Scan Center Manager

MANAGE PATIENTS

Search all columns...

+ Add Patients

Name	Patient ID	Mail ID	Mobile Number	Status	Manage Profile
Sri shaa	P100000	rohiteagle007@gmail.com	9876543214	Active	Edit
Patient #1	P100001	varunhoskere3@gmail.com	7206810456	Active	Edit
Test Ramanan	P100002	yejal49644@iamtile.com	9879799799	Active	Edit
WG	P100004	wellthgreen@gmail.com	3333333333	Active	Edit
MS	P100005	drmy3shankarmd@gmail.com	2343252352	Active	Edit

Rows per page 5

<<

< Previous

1

Next >

>>

Page 1 of 12

Copyright © Wellthgreen

The centre manager can add and manage patients. The admin can add patient by click on “Add Patient”.

Dashboard

Administration

Patient Queue

Analytics

M

MSC Manager

CMSCM2

Scan Center Manager

ADD PATIENT

Add Photo

Personal Details

Full Name *

Enter Full Name

Contact Number *

USA (+1)

Enter Phone Number

Gender *

Select Gender

Date Of Appointment *

Select a date

Patient ID *

Enter Patient ID

E-Mail *

m@example.com

Date Of Birth *


Select a date

< Back

Copyright © Wellthgreen

Submit >


The Center admin can add Patients details.



The screenshot shows the 'MANAGE PATIENTS' section of the 'easeQT' application. At the top, there are navigation tabs: Dashboard, Administration, Patient Queue, and Analytics. The 'Administration' tab is active. Below the tabs is a search bar and an 'Add Patients' button. The main area contains a table with the following columns: Name, Patient ID, Mail ID, Mobile Number, Status, and Manage Profile. The table lists five patients, all with a status of 'Active'. The 'Edit' button in the 'Manage Profile' column for the first patient, 'Sri shaa', is highlighted with a red box.

Name	Patient ID	Mail ID	Mobile Number	Status	Manage Profile
Sri shaa	P100000	rohiteagle007@gmail.com	9876543214	Active	Edit
Patient #1	P100001	varunhoskere3@gmail.com	7206810456	Active	Edit
Test Ramanan	P100002	yejal49644@iamtile.com	9879799799	Active	Edit
WG	P100004	wellthgreen@gmail.com	3333333333	Active	Edit
MS	P100005	drmy3shankarnd@gmail.com	2343252352	Active	Edit

At the bottom of the table, there is a 'Rows per page' dropdown set to 5, and pagination controls showing 'Page 1 of 12'.



The screenshot shows the 'Edit Performing Provider' form. The form is titled 'Edit Performing Provider' and has a close button in the top right corner. It contains a profile picture placeholder and the following fields:

- Patient ID:** P100000
- Status:** Active (dropdown menu)
- Personal Details:**
 - Full Name:** Sri shaa
 - E-Mail:** rohiteagle007@gmail.com
 - Contact Number:** 9876543214
 - Date Of Birth:** Select a date (calendar icon)
 - Gender:** Select Gender (dropdown menu)

At the bottom of the form, there is a 'Save Changes' button.

If the centre admin wants to change the information in a particular column, they can click the “Edit” option, make the necessary changes, and submit.

Personal Details

Full Name * E-Mail *

Contact Number * Date Of Birth *

Gender *

Appointment

Date Of Appointment *

User ID	Data	Mail
No Reports		

Copyright © Wellthgreen

The centre manager will select the date and click **Create Appointment** to book it.

Send Mail Option

Email Send Mail Option

Mail sent by admin@wellthgreen.com

Welcome, Sri shaa!

You have successfully been onboarded as a Patient at **Wellthgreen Theranostics.**

Your login credentials are as follows:

Patient ID: P100000

Email: rohteagle007@gmail.com

Password: Kindly use your password

Copyright © Wellthgreen

Send Mail Option

Email Send Mail Option

Subject

Subject

Hi Sri shaa,

Thank you for choosing MSC Scan Centre for your Breast QT (Quantitative Transmission) Imaging, also known as Breast Acoustic CT. We're looking forward to supporting you with a smooth and comfortable experience on your appointment date Sep 24, 2025

- Log in: <https://easeqremote.brightoncloudtech.com>
- Enter details: Once logged in with your email ID **rohteagle007@gmail.com** using your password.

If you have any questions or need assistance, please don't hesitate to reach out to us. We're here to make this process as seamless as possible for you. Warm regards, MSC Scan Centre Team - Contact Info

Copyright © Wellthgreen

After clicking **Create Appointment**, it redirects to the mail-sending page. There are two options:

1. Send the email directly through our application.
2. Copy the email content and send it manually through the scan centre's own email account.

PATIENT QUEUE

easeQT

easy scan transformation

Dashboard

Administration

Patient Queue

Analytics

MSC Manager

CMSCM2

Scan Center Manager

Download Dicom

Download Final Reports

Download Patient Consent

Clear All Filters

Search all columns...

	Date	Scan Centre	Patient Name	Consent	Form	Tech Form	DICOM		DICOM	Old Report	QT Report	Report Status	Report Delivery	Assign	Pending Remarks	Remarks
							L	R								
<input type="checkbox"/>	Sep 17, 2025	MSC	Testing12	View	Dc - View	View	(1)		View DICOM		View	SIGNED OFF (A)	-		t dicom pen...	<div><div></div><div></div></div>
<input type="checkbox"/>	Sep 18, 2025	MSC	Testing12	View	S - View	View	(1)		View DICOM		-	REVIEWED 2	-	Assign	R1 R2 stuck	<div><div></div><div></div></div>
<input type="checkbox"/>	Sep 22, 2025	MSC	Testing12	View	S - View	View	(1)		View DICOM		View	SIGNED OFF (A)	-	Assign	No remarks...	<div><div></div><div></div></div>
<input type="checkbox"/>	Sep 18, 2025	MSC	sofia	View	S - View	View	(1)		View DICOM		-	DRAFT	-	Assign	No remarks...	<div><div></div><div></div></div>
<input type="checkbox"/>	Sep 19, 2025	MSC	sofia	View	S - View	View	(1)		View DICOM		-	PREDRAFT	-	Assign	No remarks...	<div><div></div><div></div></div>
<input type="checkbox"/>	Sep 19, 2025	MSC	U S	View	S - View	View	(1)		View DICOM		-	REVIEWED 2	-	Assign	No remarks...	<div><div></div><div></div></div>

Rows per page

50

<<

< Previous

1


Next >

>>

Page 1 of 1

Copyright © Wellthreen

In the patient queue, the scan centre admin can view and download patient reports. They can also edit pending remarks. They can upload and download patient's old reports.



Dashboard

Administration

Patient Queue

Analytics

t

testing scan centre manager

CTSCM1

Scan Center Manager

Download Dicom

Download Final Reports

Download Patient Consent

Clear All Filters

Search all columns...

	Date	Scan Centre	Patient Name	Consent	Form	Technician Form	DICOM		Report	Report Status	Report Delivery	Assigned	Pending Remarks	Remarks
							L	R						
<input type="checkbox"/>	8/7/2025	TSC	testing patient	View	Dc - View	View	(1)	(1)	-	YET TO REPORT	-	Assign	No remarks yet	+
<input type="checkbox"/>	8/6/2025	TSC	testing patient	View	S - View	View	(1)	(1)	View	SIGNED OFF	Send Mail	WGPP2	No remarks yet	+
<input type="checkbox"/>	8/8/2025	TSC	sofia	View	S - View	View	(2)	-	-	YET TO REPORT	-	Assign	sdddsds	+
<input type="checkbox"/>	8/9/2025	TSC	sofia	View	S - View	View	(2)	(1)	View	SIGNED OFF	Send Mail	Assign	No remarks yet	+

Rows per page 50

<<

< Previous

1

Next >

>>

Page 1 of 1

Dashboard

Administration

Patient Queue

Analytics

?

M

MSC Manager

CMSCM2

Scan Center Manager

Download Dicom

Download Final Reports

Download Patient Consent

Clear All Filters

Search all columns...

	Date	Scan Centre	Patient Name	Consent	Form	Tech Form	DICOM	DICOM	Old Report	QT Report	Report Status	Report Delivery	Assign	Pending Remarks	Remarks
	Sep 17, 2025	MSC	Testing12	View	Dc - View	View	(1)	-	View DICOM		View	SIGNED OFF (A)	Send Mail		I dicom pen...
	Sep 18, 2025	MSC	Testing12	View	S - View	View	(1)	-	View DICOM		View	REVIEWED 2		Assign	R1 R2 stuck
	Sep 22, 2025	MSC	Testing12	View	S - View	View	(1)	-	View DICOM		View	SIGNED OFF (A)	Send Mail	Assign	No remarks...
	Sep 18, 2025	MSC	sofia	View	S - View	View	(1)	-	View DICOM		View	DRAFT		Assign	No remarks...
	Sep 19, 2025	MSC	sofia	View	S - View	View	(1)	-	View DICOM		View	PREDRAFT		Assign	No remarks...
	Sep 19,								View			REVIEWED			

Rows per page

50

<<

< Previous

1

Next >

>>

Page 1 of 1

Copyright © Wellthreen

Dashboard

Administration

Patient Queue

Analytics

t

testing scan centre manager

CTSCM1

Scan Center Manager

Download Dicom

Download Final Reports

Download Patient Consent

Clear All Filters

Search all columns...

	Date	Scan Centre	Patient Name	Consent	Form	Technician Form	DICOM	DICOM	Report	Report Status	Report Delivery	Assigned	Pending Remarks	Remarks
	8/7/2025	TSC	testing patient	View	Dc - View							Assign	No remarks yet	
	8/6/2025	TSC	testing patient	View	S - View						Send Mail	WGPP2	No remarks yet	
	8/8/2025	TSC	sofia	View	S - View							Assign	sdddsds	
	8/9/2025	TSC	sofia	View	S - View						Send Mail	Assign	No remarks yet	

Rows per page

50

<<

< Previous

1

Next >

>>

Page 1 of 1

Copyright © Wellthreen

Select Email Recipients

Select recipients to proceed with sending the email.

Choose recipient

Send

Dashboard

Administration

Patient Queue

Analytics

t

testing scan centre manager

CTSCM1

Scan Center Manager

Download Dicom

Download Final Reports

Download Patient Consent

Clear All Filters

Search all columns...

	Date	Scan Centre	Patient Name	Consent	Form	Technician Form	DICOM	DICOM	Report	Report Status	Report Delivery	Assigned	Pending Remarks	Remarks
	8/7/2025	TSC	testing patient	View	Dc - View							Assign	No remarks yet	
	8/6/2025	TSC	testing patient	View	S - View						Send Mail	WGPP2	No remarks yet	
	8/8/2025	TSC	sofia	View	S - View							Assign	sdddsds	
	8/9/2025	TSC	sofia	View	S - View						Send Mail	Assign	No remarks yet	

Rows per page

50

<<

< Previous

1

Next >

>>

Page 1 of 1

Copyright © Wellthreen

Select Email Recipients

Select recipients to proceed with sending the email.

Choose recipient


Patient

Scan Center Manager

Patient and Scan Center Manager

After the doctor completes the report, the centre manager can send the final report to the patient.

ANALYTICS



In Analytics, the centre manager can view the performance and activity analytics of all employees within their scan centre.

GENERAL TIPS

1. If you open something in edit mode, save it before exiting-otherwise, others won't be able to use it.
2. After addressing pending remarks, leave another remark saying "Done" and set Assign to as None.
3. The Filter option is available for sorting.
4. For bulk downloads, use the checkboxes to select patients and then download in bulk.
5. To add a pending remark, click the + icon.
6. When viewing a report, always click Back to exit-otherwise, the case you opened may freeze.
7. Use Image Preview to add images to a report.
8. To refresh the page, press Ctrl + Shift + R.
9. For DICOM downloads and voice dictation, use Chrome or Edge. (For Firefox, click Save. Do not use Internet Explorer.)
10. Always log out properly-do not just close the browser or tab.