




ease @T


SOFT, SAFE, SYNCED BREAST HEALTH

SCRIBE LOGIN USER MANUAL



Once the technician completes the form and uploads the DICOM, the scribe will begin preparing the report.

The scribe signs in using the email and password provided by WellthGreen.



A verification code will be generated and sent for verification.

A screenshot of a web browser displaying the 'EaseQT Platform' User Consent Form. The form is titled 'User Consent Form' and 'EaseQT Platform'. It contains four sections: 1. Consent to Platform Use, 2. Data Privacy & Security, 3. Confidentiality Agreement, and 4. Data Consent. Each section contains a paragraph of text explaining the terms of use and consent. The background of the browser window shows a collage of medical images and a person using a computer.

easeQT
SOFT, SAFE, SYNCED BREAST HEALTH

User Consent Form

EaseQT Platform


1. Consent to Platform Use
I hereby acknowledge that I am voluntarily registering to use the EaseQT platform for purposes directly related to clinical imaging, workflow management, or reporting.

2. Data Privacy & Security
I understand that my personal information and any patient-related data accessed through this platform may be securely stored and processed in compliance with global privacy laws (including HIPAA, DPA 2021, and GDPR where applicable). I agree not to share, download, or misuse any protected information obtained via the platform.

3. Confidentiality Agreement
I agree to maintain full confidentiality of any data I access, including but not limited to patient records, reports, clinical communications, and platform operations. Any breach of confidentiality may result in legal and professional consequences.

4. Data Consent
I consent to the transfer, storage, and processing of data as necessary for platform function and continuity of care. I understand that all such data will be handled securely and only by authorized personnel.

The scribe will sign the consent form.

A screenshot of the EaseQT Password Reset page. The page features the EaseQT logo and the tagline 'SOFT, SAFE, SYNCED BREAST HEALTH'. Below the logo, it says 'Reset Your Password' and 'E-Mail: techniciancentre@gmail.com'. There are two password input fields: 'New Password' and 'Confirm Password'. The 'New Password' field has a dropdown menu with options: 'At least 8 characters', 'Contains a number', and 'Special character (!@#%&*^&*)'. The 'Confirm Password' field has a dropdown menu with the option 'Confirm password'. A 'Reset Password' button is located at the bottom of the form. The background of the page shows a collage of medical images and a person using a computer.

easeQT
SOFT, SAFE, SYNCED BREAST HEALTH

Reset Your Password

E-Mail: techniciancentre@gmail.com

New Password

Password

At least 8 characters
Contains a number
Special character (!@#%&*^&*)


Confirm Password

Confirm password

Reset Password


The scribe needs to reset their password before logging in.

DASHBOARD




In the dashboard, the scribe can view the patient brochure, technician guidelines, consent form and invoices.

1. Invoices




Clicking on the 'Admin' button triggers a popup window.



Clicking the 'Invoices' button open to the Invoice page.



This section displays a list of invoices generated by you or the admin.



If you need an invoice, select the desired month from the enabled options and click the 'Create Report' button to proceed to the invoice creation page.

The marked fields are auto-populated with data from the system and cannot be edited. Unmarked fields are editable. Once all the required details are filled in, click the 'Generate Report' button to create the invoice.

In the patient queue, the scribe can view the patient list, patient forms, technician view, and DICOM files. They can fill out reports, check the report status, assign the case to another user and add pending remarks if needed. And they can multi select and download DICOM. They can upload and download patient's old reports.

Time stamp

Created by	Start	End	Status
sofia	18 Sept 2025 16:45	18 Sept 2025 04:16	Intake Form Fill Complete
Technician	18 Sept 2025 04:17	18 Sept 2025 04:18	Technology Form Fill
Wellthgreen Scribe	18 Sept 2025 04:18		

Status saver

Load Template Reset to Default

Predraft Draft

Reviewed 1 Correct Reviewed 1 Edit

Reviewed 2 Correct Reviewed 2 Edit

Sign Off

A. PATIENT HISTORY

Report Preview

Normal B I U [List Icon] [Align Icon] [Indent Icon] [Decrease Indent Icon] A [Link Icon] [Image Icon]

Indication: Screening : Routine Annual Checkup.
History:
Reason for having this QT scan (Patient Stated): sss.
Patient is asymptomatic.
Imaging:
• sss.

B. BREAST IMPLANT DETAILS

Breast Implants ☐ Present ☒ Absent

Report Preview

Copyright © Wellthreen

In the report section, the scribe can review and modify reports completed by the technician. A "Pre-draft" button is available for the scribe's use. Once the technician has completed the form and uploaded the DICOM files, the scribe can either select predefined options or directly edit the report using the text editor.

Other reports

S.No	Report Title	File
No Reports Found		

Other Old reports

S.No	Report Title	File
1	Thermogram	View

Final Report Preview

Ease QT 10.10 Auto Report ☒


Normal B I U [List Icon] [Align Icon] [Indent Icon] [Decrease Indent Icon] A [Link Icon] [Image Icon]

NAME	sofia	DOB	August 4, 1997 (28 y)
GENDER	F	SCAN CENTER	MSC, Bangalore
USERID	P100016	DATE OF VISIT	Sep 18, 2025


QT ULTRASOUND BREAST IMAGING
Indication: Screening : Routine Annual Checkup.
History:
Reason for having this QT scan (Patient Stated): sss.
Patient is asymptomatic.
Imaging:
• sss.

TECHNIQUE: Transmission and reflection multiplanar 3-dimensional ultrasound imaging of both breasts was performed using the QT Ultrasound Series 2000 Model-A scanner.

Copyright © Wellthreen



The users can see uploaded Patient's old reports and their current report.



The scribe can view the full report here. If they want to make changes or load a template, they need to click “Unsync” to switch it to “Sync,” make the necessary changes, and update the report.

ANALYTICS



In Analytics, the scribe can view their case count, the case count for the previous six months.

GENERAL TIPS

1. If you open something in edit mode, save it before exiting-otherwise, others won't be able to use it.
2. After addressing pending remarks, leave another remark saying "Done" and set Assign to as None.
3. The Filter option is available for sorting.
4. For bulk downloads, use the checkboxes to select patients and then download in bulk.
5. To add a pending remark, click the + icon.
6. When viewing a report, always click Back to exit-otherwise, the case you opened may freeze.
7. Use Image Preview to add images to a report.
8. To refresh the page, press Ctrl + Shift + R.
9. For DICOM downloads and voice dictation, use Chrome or Edge. (For Firefox, click Save. Do not use Internet Explorer.)
10. Always log out properly-do not just close the browser or tab.