



ease @T

SOFT, SAFE, SYNCED BREAST HEALTH


WELLGREEN MANAGER LOGIN USER MANUAL



The wellthgreen manager will sign in using email and password provided by wellgreen.




The passcode will be generated for verification.



The screenshot shows a web browser window with the title "User Consent Form" and the subtitle "EaseQT Platform". The form is divided into four sections:

- 1. Consent to Platform Use**
I hereby acknowledge that I am voluntarily registering to use the EaseQT platform for purposes directly related to clinical imaging, workflow management, or reporting.
- 2. Data Privacy & Security**
I understand that my personal information and any patient-related data accessed through this platform may be securely stored and processed in compliance with global privacy laws (including HIPAA, DPA 2021, and GDPR where applicable). I agree not to share, download, or misuse any protected information obtained via the platform.
- 3. Confidentiality Agreement**
I agree to maintain full confidentiality of any data I access, including but not limited to patient records, reports, clinical communications, and platform operations. Any breach of confidentiality may result in legal and professional consequences.
- 4. Data Consent**
I consent to the transfer, storage, and processing of data as necessary for platform function and continuity of care. I understand that all such data will be handled securely and only by authorized personnel.


The wellthgreen manager needs to sign the form.




The screenshot shows a web browser window with the title "Reset Your Password" and the subtitle "EaseQT". The page features the EaseQT logo at the top, which includes the text "SOFT, SAFE, SYNCHRONIZED BREAST HEALTH". Below the logo, the text "Reset Your Password" is displayed, followed by the email address "E-Mail: gokulhk278@gmail.com". The form includes two password input fields: "New Password" and "Confirm Password", both with a strength indicator. Below the "New Password" field, there are three checkboxes indicating password requirements: "At least 8 characters", "Contains a number", and "Special character (!@#%&*)". At the bottom of the form, there is a "Reset Password" button.

The wellthgreen manager needs to reset their password.

DASHBOARD



In the dashboard, the wellthgreen manager can view the patient brochure, general guidelines, disclaimer, patient consent form, technician guidelines, and technician consent form. They can also view the radiologist training material, radiologist consent form, and invoices.



The screenshot shows a web form titled "Invoice" with the easeQT logo at the top left. The form is divided into three main sections: "From", "Billing Period", and "Invoice Details".

- From:** Includes input fields for "Wellth Green Health Care Pvt. Ltd.", "Phone Number", "Email", "PAN", "GST", and "Address".
- Billing Period:** Includes date pickers for "From" (June 1st, 2025) and "To" (June 30th, 2025), a dropdown for "UPI", and an input field for "UPI ID".
- Invoice Details:** Includes a table with "Quantity" (0) and "Amount (INR)" (7000), a "Total" field showing "0/-", and a "Generate Report" button.

A "Back" button is located at the top left of the form area.

The screenshot shows a modal window titled "Invoice" with the easeQT logo at the top left. The modal is overlaid on a dashboard background. The modal contains the following elements:

- Buttons:** "Scan Center", "User", and "Edit Amount".
- Form:** A dropdown menu labeled "Choose Scan Center" with a checkmark icon.
- Message:** "No Invoice Found" displayed on the right side of the modal.

The background dashboard shows a "Welcome Admin" message, a user profile for "Admin AD100001", and various system icons.

For invoice generation, the admin can select either a scan centre or a user, choose the desired month, and click “Create Report.” A form will appear where the required details need to be entered, after which the report will be generated.

Invoice

← Back

From

Wellthgreen Performing Provider

9842653436

wgpp@gmail.com

PAN

GST

Address

Billing Period

June 1st, 2025

To

June 30th, 2025

UPI

UPI ID

Invoice Details

Quantity	Amount (INR)
0	3000678
Total 0/-	

To

Wellth Green Health Care Pvt. Ltd.

Address

Generate Report

The marked fields are auto-populated with data from the system and cannot be edited. Unmarked fields are editable. Once all the required details are filled in, click the 'Generate Report' button to create the invoice.

Invoice

Scan Center User Edit Amount


Scan Center Amount (INR) User Amount (INR)

677 3000678


Save

The manager can fix the amount for scan centre and user.

ADMINISTRATION



In Administration, the wellthgreen manager can add and manage each user under WellthGreen Administration.



The wellthgreen manager can add performing provider by click on 'Add performing provider'.

Dashboard
Administration
Patient Queue
Analytics

wellth green manager
WGM1

ADD WELLTHGREEN PERFORMING PROVIDER

Add Photo

Personal Details

Full Name *

E-Mail *

Aadhar *

Driver's License *

Contact Number *
IN (+91)

Date Of Birth *

PAN *

Professional Details

MBBS Registration Number *

Specialization *

MD Registration Number *

Upload CV *

Back
Submit

The wellthgreen manager can add the personal and professional details required to add a performing provider. They can also attach the person's Aadhaar, PAN, driver's license, MBBS registration number, MD registration etc.

Dashboard
Administration
Patient Queue
Analytics

wellth green manager
WGM1

MANAGE WELLTHGREEN PERFORMING PROVIDER

+ Add Performing Provider

Name	WGPP ID	Mail ID	Mobile Number	Status	Manage Profile
scan performing provider	WGPP6	scanprovider@gmail.com	1234567899	Active	Edit
scan performing provider	WGPP5	scanpp@gmail.com	1234567444	Active	Edit
S Performing Provider	WGPP4	sperformingprovider@gmail.com	1234567890	Active	Edit
Rohith	WGPP3	samplewgpp2@gmail.com	9123456789	Active	Edit
sampleWGPP1	WGPP2	samplewgpp1@gmail.com	9842653420	Active	Edit

Rows per page: 5

Previous
1
Next

Page 1 of 2

The screenshot shows a web application interface for managing radiologists. A modal window titled "Edit Radiologist" is open, displaying a form for updating a radiologist's information. The form includes a profile picture placeholder, a "Radiologist ID" field with the value "WGPP6", and a "Status" dropdown menu currently set to "Active". Below these is a section titled "Personal Details" which contains several fields: "Full Name" (with the value "scan performing provider"), "E-Mail" (with the value "scanprovider@gmail.com"), "Aadhar" (with an "Upload Aadhar" button), "Driving License" (with an "Upload Driving License" button), "Contact Number" (with the value "1234567899"), "Date Of Birth" (with the value "July 1st, 2025" and a calendar icon), and "PAN" (with an "Upload PAN" button). The background shows a sidebar with a "MANAGE" section and a table of radiologists.

Edit Radiologist

Radiologist ID: WGPP6

Status: Active

Personal Details

Full Name: scan performing provider

E-Mail: scanprovider@gmail.com

Aadhar: Upload Aadhar

Driving License: Upload Driving License

Contact Number: 1234567899

Date Of Birth: July 1st, 2025

PAN: Upload PAN

If the wellthgreen manager wants to change the information in a particular column, they can click the “Edit” option, make the necessary changes, and submit. If they want to deactivate a user, they can change the user’s status from “Active” to “Inactive.”

The same process applies for radiologists, and scribes.

PATIENT QUEUE

	Date	Scan Centre	Patient ID	Patient Name	Consent	Form	Tech Form	DICOM	Old Report	QT Report	Report Status	Assign	Pending Remarks	Remarks
<input type="checkbox"/>	Sep 17, 2025	MSC	P100006	Testing12	View	Dc - View	View	(1) -	View DICOM		View	SIGNED OFF (A)	AD1000C	No remarks...
<input type="checkbox"/>	Sep 18, 2025	MSC	P100006	Testing12	View	S - View	View	(1) -	View DICOM		-	REVIEWED 2	Assign	R1 R2 stuck
<input type="checkbox"/>	Sep 22, 2025	MSC	P100006	Testing12	View	S - View	View	(1) -	View DICOM		View	SIGNED OFF (A)	Assign	No remarks...
<input type="checkbox"/>	Sep 18, 2025	MSC	P100016	sofia	View	S - View	View	(1) -	View DICOM		-	DRAFT	Assign	No remarks...

Rows per page: 50

« < Previous 1 Next > »

Page 1 of 1

Copyright © Wellthgreen

← Back

Patient Name: Testing12
Patient ID: P100006
Gender: Female
DOB: January 4, 1972 (53 Y)
Patient Form: Dc Form
Date: 17 Sept 2025
Scan Center: MSC

Other reports

S.No	Report Title	File
No Reports Found		

Other Old reports

S.No	Report Title	File
1	Thermogram	View

Time stamp

Created by	Start	End	Status
Invalid	17	Patient	

Final Report Preview

Ease QT 10:10 Auto Report

Normal

B

I

U

GENDER	F	SCAN CENTER	MSC, Bangalore
USERID	P100006	DATE OF VISIT	Sep 17, 2025

QT ULTRASOUND BREAST IMAGING

Indication: Comparative Study as follow up to a prior QT scan.

History:

Reason for having this QT scan (Patient Stated): 545.

Breast implant present on the left, silicone type, since 12months, right, silicone type, since 12months.

Current breast symptoms:

Breast pain in the right breast 4 o'clock since , so


TECHNIQUE: Transmission and reflection multiplanar 3-dimensional ultrasound imaging of both breasts was performed using the QT Ultrasound Series 2000 Model-A scanner. Breast density was determined using the Quantitative Breast Density calculator. Images were reviewed in the QTviewer v2.6.2. The nipple-areolar complex, skin, Cooper's ligaments, breast fat distribution, penetrating arteries and veins, glandular and ductal tissues were evaluated. Images were reviewed in coronal, transaxial and sagittal planes.

BREAST IMPLANTS:

The QT scan shows bilateral similar implants which are subpectoral in position, with

In patient queue, the wellthgreen manager can view the patient list, patient forms, technician view, and DICOM files. the wellthgreen manager can view the report, check the report status, assign the case to another user, and add pending remarks if needed and add pending remark and can view previous remark also. They can multi select and can download the DICOM and Final reports. The manager can view the report.

ANALYTICS



In analytics, The wellthgreen manager can view all users and scan centres analytics.