



ease @T



**SOFT, SAFE, SYNCED BREAST HEALTH**

## SCRIBE LOGIN USER MANUAL



Once the technician completes the form and uploads the DICOM, the scribe will begin preparing the report.

The scribe signs in using the email and password provided by WellthGreen.



A passcode will be generated and sent for verification.



The screenshot shows a web browser window with the title "User Consent Form" and the subtitle "EaseQT Platform". The form is divided into four sections, each with a heading and a paragraph of text. The background of the browser window shows a collage of medical images and a person in a lab coat.

**easeQT**  
SOFT, SAFE, SYNCED BREAST HEALTH

### 1. Consent to Platform Use

I hereby acknowledge that I am voluntarily registering to use the EaseQT platform for purposes directly related to clinical imaging, workflow management, or reporting.

### 2. Data Privacy & Security

I understand that my personal information and any patient-related data accessed through this platform may be securely stored and processed in compliance with global privacy laws (including HIPAA, DPA 2021, and GDPR where applicable). I agree not to share, download, or misuse any protected information obtained via the platform.


### 3. Confidentiality Agreement

I agree to maintain full confidentiality of any data I access, including but not limited to patient records, reports, clinical communications, and platform operations. Any breach of confidentiality may result in legal and professional consequences.

### 4. Data Consent

I consent to the transfer, storage, and processing of data as necessary for platform function and continuity of care. I understand that all such data will be handled securely and only by authorized personnel.

The scribe will sign the consent form.



The screenshot shows a web browser window with the title "Reset Your Password" and the subtitle "EaseQT Platform". The form is divided into two main sections: a left sidebar with a collage of medical images and a right main content area with a password reset form. The background of the browser window shows a collage of medical images and a person in a lab coat.

**easeQT**  
SOFT, SAFE, SYNCED BREAST HEALTH

### Reset Your Password

E-Mail: [techniciancentre@gmail.com](mailto:techniciancentre@gmail.com)

New Password

Password

At least 8 characters  
Contains a number  
Special character (!@#%&\*)


Confirm Password

Confirm password

**Reset Password**


The scribe needs to reset their password before logging in.

## DASHBOARD




In the dashboard, the scribe can view the patient brochure, technician guidelines, consent form and invoices.

### 1. Invoices




Clicking on the 'Admin' button triggers a popup window.



Clicking the 'Invoices' button open to the Invoice page.




This section displays a list of invoices generated by you or the admin.




If you need an invoice, select the desired month from the enabled options and click the 'Create Report' button to proceed to the invoice creation page.

The marked fields are auto-populated with data from the system and cannot be edited. Unmarked fields are editable. Once all the required details are filled in, click the 'Generate Report' button to create the invoice.



Generated invoices are listed here. If you need a copy, click the 'Download' button to download your invoice.



## Invoice Preview

## PATIENT QUEUE

easeQT

easy soft. intelligent health

Dashboard

Patient Queue

Analytics

0

W

Wellthgreen Scribe

WGS4

WellthGreen Scribe

Download Dicom

Clear All Filters

Search all columns...

<input type="checkbox"/>	Date	Scan Centre	Patient ID	Patient Name	Form	Tech Form	DICOM		DICOM	Old Report	QT Report	Report Status	Assign	Pending Remarks	Remarks
							L	R							
<input type="checkbox"/>	Sep 17, 2025	MSC	P100006	Testing12	Dc - View	View	(1)	-	View DICOM		View	SIGNED OFF (A)	<div><div></div></div>	No remarks...	<div><div>+</div><div></div></div>
<input type="checkbox"/>	Sep 18, 2025	MSC	P100006	Testing12	S - View	View	(1)	-	View DICOM		View	REVIEWED 2	<div>Assign <div></div></div>	R1 R2 stuck	<div><div>+</div><div></div></div>
<input type="checkbox"/>	Sep 22, 2025	MSC	P100006	Testing12	S - View	View	(1)	-	View DICOM		View	SIGNED OFF (A)	<div>Assign <div></div></div>	No remarks...	<div><div>+</div><div></div></div>

Rows per page

50

<<

< Previous

1

Next >

>>

Page 1 of 1

Copyright © Wellthgreen

In the patient queue, the scribe can view the patient list, patient forms, technician view, and DICOM files. They can fill out reports, check the report status, assign the case to another user, and add pending remarks if needed. And they can multi select and download DICOM. They can upload and download patient's old reports.

Time stamp

Created by	Start	End	Status
sofia	2025-09-16 16:45	2025-09-16 04:16	Completed
Technician	2025-09-18 04:17	2025-09-18 04:18	Technology Form Fill
Wellthgreen Scribe	2025-09-18 04:18	-	-

Status saver

Load Template

Reset to Default

Predraft

Draft

Reviewed 1 Correct

Reviewed 1 Edit

Reviewed 2 Correct

Reviewed 2 Edit

Sign Off

A. PATIENT HISTORY

Report Preview

Normal

Indication: Screening : Routine Annual Checkup.

History:

Reason for having this QT scan (Patient Stated): sss.

Patient is asymptomatic.

Imaging:

• sss.

B. BREAST IMPLANT DETAILS

Breast Implants

☐ Present

☒ Absent

Report Preview

In the report section, the scribe can review and modify reports completed by the technician. A "Pre-draft" button is available for the scribe's use. Once the technician has completed the form and uploaded the DICOM files, the scribe can either select predefined options or directly edit the report using the text editor.





**Final Report Preview** Ease QT 10.10 Auto Report ☐

**Time stamp**

Created by	Start	End	Status
sofia	2025-09-18 16:45	2025-09-18 04:16	Fill
Technician	2025-09-18 04:17	2025-09-18 04:18	Completed
Wellthgreen Scribe	2025-09-18 04:18	-	-

**Status saver**

Load Template

**Patient Information:**

NAME	sofia	DOB	August 4, 1997 (28 y)
GENDER	F	SCAN CENTER	MSC, Bangalore
USERID	P100016	DATE OF VISIT	Sep 18, 2025

**QT ULTRASOUND BREAST IMAGING**  
 Indication: Screening : Routine Annual Checkup.  
 History:  
 Reason for having this QT scan (Patient Stated): sss.  
 Patient is asymptomatic.  
 Imaging:  
 • sss.

**TECHNIQUE:** Transmission and reflection multiplanar 3-dimensional ultrasound imaging

Copyright © Wellthgreen

The scribe can view the full report here. If they want to make changes or load a template, they need to click “Unsync” to switch it to “Sync,” make the necessary changes, and update the report.

## ANALYTICS



In Analytics, the scribe can view their case count, the case count for the previous six months.